
Licel Transient Recorder and Ethernet-Controller

Programming Manual

Licel GmbH

March 31, 2023

Contents

1	Introduction	5
2	Installation	6
2.1	Ethernet software	6
2.1.1	Required software	6
2.2	Software installation	6
2.2.1	LabVIEW Sources	6
2.2.2	C-sources	6
3	Program Organization	6
4	Modules	6
4.1	licel_tcpip	6
4.2	licel_tr/licel_tr_tcpip	7
4.2.1	Controller	7
4.3	licel_data	13
4.4	licel_util	13
4.5	APD functions	14
4.6	PMT functions	15
4.7	Timing functions	15
4.8	Security functions	16
4.9	Power Meter	16
4.10	Bore Alignment	16
4.11	Wind Acquisition	17
4.12	Function arguments	18
4.13	Timing Parameter explanation	22
5	Transient recorder	24
5.1	Low Level Commands	24
5.2	Memory organization	25
5.3	Raw Data to Physical Value Conversion	27
5.4	Raw Squared data to signal standard deviation	28
5.5	Acquisition Low Level Description	29
5.6	High Performant Read Operations	30
6	Network security	31
7	C - Example Programs	32
7.1	The script	33
7.2	start	33
7.3	shot	33
7.4	readout	33
7.5	Example Program Configuration File	33
7.6	File format	34
7.7	SampleAcquis for Ethernet Applications	36
7.8	MPUSH Mode Sample Acquisition	36
7.9	APD Control example	38
7.10	Wind Sample Acquisition	38
7.11	Network management utilities	40
7.11.1	Getting Started	40
7.11.2	Set Fixed IP Address	40
7.11.3	Activate DHCP Mode	40
7.11.4	SecureModeEnable	40
7.11.5	SecureModeDisable	40

8 Appendix VB6/VB.net Programming	40
8.1 Sample applications	40
8.1.1 Control Overview	40
8.1.2 MultipleChannel	40
8.1.3 PushModeDemo	41
9 Appendix - Obsolete DIO32HS references	41
9.1 Operation Principles	41
9.2 Hardware Requirements	41
9.2.1 Further References	41
9.3 Windows - DIO-32HS	41
9.3.1 Required software	41
9.3.2 NI-DAQ-Setup	42
9.3.3 Interface card installation	42
9.3.4 Verification	43
9.4 Linux - PCI-DIO-32HS	43
9.4.1 Kernel preparation	43
9.4.2 Necessary Files	43
9.4.3 Driver tests and card installation	43
9.4.4 Changes to /etc/modules.conf	44
9.4.5 Changes to /etc/rc.d/rc.local	44
9.4.6 Directory structure	44
9.4.7 LabView + comedi	44
9.4.8 licel_nidaq/licel_re	44

1 Introduction

Currently the Transient Recorder can be controlled by two ways:

- Via a Ethernet interface. The Ethernet interface has then its own parallel bus interface to the Transient Recorders.
- Via the parallel interface cards from National Instruments, this option is deprecated. It only works till Win-XP. There is no driver support under NIDAQmx and this excludes Windows Vista Windows 7 as operating systems. see the Appendix for historical references if you encounter such a system.

The software to control the Licel Transient Recorder will

- Control one or more Transient Recorders
- Ensure software portability between different operating systems
- Readout the Transient Recorders at high data transfer rates.

The software for the Ethernet Controller will in addition control

- The APD module
- The PMT module
- The Trigger module

The following target systems are supported:

- LabVIEW from National Instruments
 - WinXX, Linux
 - Mac-OS
- MS Visual C and gcc for the Ethernet Controller software.
- A Visual Basic/VB.net module for the Ethernet Controller software. This module is sold separately.
- *obsolete* NIDAQ-Library from National Instruments
 - WinXX with MS Visual C
 - WinXX with MS Visual Basic
- *obsolete* A COMEDI based Linux driver for gcc

The software is able to:

- configure the Transient Recorders
- start the data acquisition
- stop the data acquisition
- query the Transient Recorder status
- readout the Transient Recorder
- convert the binary data to quantities with physical units.

The Ethernet based software additionally is able to:

- Set the PMT high voltage and read the PMT status.
- Set the APD high voltage, Activate the APD thermo electrical cooler and read the APD Status back.
- Set the delays on the trigger generator and to activate separately the trigger lines.

2 Installation

2.1 Ethernet software

2.1.1 Required software

In order to use the C,Basic or LabView routines under Windows or Linux, one needs a working installation of

- Visual C++ 6.0 or gcc(Linux)
- Visual Basic 6.0 .
- LabView 10.x or higher, the Wind-software requires LabVIEW 2018 or higher

The controlling PC should have a working network connection. The details of the of setting up the Ethernet Controller are outlined in the [Ethernet Manual - Chapter 4](#).

2.2 Software installation

2.2.1 LabVIEW Sources

please follow the instructions in the [Ethernet Manual - Chapter 3](#)


2.2.2 C-sources

Please download the sources at http://licel.com/download/c-files/Ethernet/licel_tcpip_C_driver.zip
And unzip the sources.

3 Program Organization

The software has layered structure where the low-level routines are encapsulated in `licel_nidaq.c` or `licel_nidaq.bas`. These routines access NI-DAQ. Under Linux the corresponding file is `licel_re.c`. This file calls the corresponding comedilib functions. Above this layer are `licel_tr.c` or `licel_tr.bas`. At this layer the functional tasks are translated into low-level commands.

For the Ethernet Controller the basic communication routines are inside `licel_tcpip.cpp` Above this level are the `licel_tr_tcpip.cpp`. At this layer the functional tasks are translated into low-level ASCII commands. Above this layer are the application functions.



4 Modules

4.1 licel_tcpip

openConnection open the connection to specified host at a specified port

```
SOCKET openConnection(const char* host, int port);
```

```
Public Function openConnection(ByVal sHost As String, ByVal iPort As Integer) As TcpClient
```

open the connection to specified host at a specified port, stays silent if it fails

```
Public Function openConnection(ByVal sHost As String, ByVal iPort As Integer, ByVal silent As Boolean) As TcpClient
```

openSecureConnection open the connection to specified host at a specified port when the access to the Controller is Limited (see [Network Security](#) for details)

```
SOCKET openSecureConnection(const char* sHost, int iPort, const char* connectionPasswd);
```

closeConnection close the specified connection

```
int closeConnection(SOCKET s);
```

```
Public Function closeConnection(ByVal client As TcpClient) As Integer
```

writeCommand write a string to the tcpconnection, append the terminating CRLF

```
int writeCommand(SOCKET s, const char* command);
```

```
Public Function writeCommand(ByVal client As TcpClient, ByVal command As String) As Integer
```

readResponse Read a ASCII response from the Controller. Except for binary data transfer the Controller response is a short string indicating whether the action could be performed or not. This response is terminated by a CRLF. This routine will read till it encounters a CRLF or if the amount of chars would exceed maxLength.

```
int readResponse(SOCKET s, char* response, int maxLength, int nTimeoutMillisec);
```

```
Public Function readResponse(ByVal client As TcpClient, ByRef response As String, ByVal maxlength As Integer, ByVal nTimeoutMillisec As Integer) As Integer
```

ReadArray Read a binary response from the Controller.

```
int ReadArray(SOCKET s, unsigned char *array, unsigned long points, int nTimeoutMillisec);
```

```
Public Function ReadArray(ByVal client As TcpClient, ByRef Data() As Byte, ByVal points As Long, ByVal nTimeoutMillisec As Integer)
```

GetPushData Read mpush data from the controller. Returns the number of bytes read in.

```
int GetPushData(SOCKET s, unsigned char* pushBuffer, int pushBuffer, int numBytes);
```

4.2 licel_tr/licel_tr_tcpip

There is a major difference in the programming model between the DIO-32 and the Ethernet version, the Ethernet version first selects the Transient Recorder and all commands (issued later) are addressed to this Transient Recorder. There are separate functions available for tasks such as starting all Transient Recorders. Before they can be called a list of Transient Recorders should be selected, the functions indicated by `Multiple` in the function name, are used for these purposes.

4.2.1 Controller

Licel_TCPIP_ActivateDHCP activate the DHCP mode on the Controller

```
int Licel_TCPIP_ActivateDHCP(SOCKET s, int iPort, const char* passwd);
```

```
Public Function Licel_TCPIP_ActivateDHCP(ByVal client As TcpClient, ByVal iPort As Integer, ByVal passwd As String) As Integer
```

Licel_TCPIP_SetIPParameter Configure the Controller for static IP configuration. Set the new IP address, the basic port number, the subnet mask and the gateway

```
int Licel_TCPIP_SetIPParameter(SOCKET s, char* newHost, char* mask,
int newPort, char* gateway, char* passwd);
```

```
Public Function Licel_TCPIP_SetIPParameter(ByVal client As TcpClient,
ByVal newHost As String, ByVal mask As String, ByVal newPort As Integer,
ByVal gateway As String, ByVal passwd As String) As Integer
```

Licel_TCPIP_GetID Get the identification string from the Controller

```
int Licel_TCPIP_GetID(SOCKET s, char* buffer, int bufferLength);
```

```
Public Function Licel_TCPIP_GetID(ByVal client As TcpClient, ByRef buffer
As String, ByVal bufferLength As Integer) As Integer
```

Licel_TCPIP_GetCapabilities Get the available subcomponents of the Controller like:

```
TR      for controlling Transient Recorder
APD     for APD remote control
PMT     for PMT remote control
TIMER   for the trigger timing Controller
CLOUD   for Transient Recorder Controller cloud mode
BORE    Boresight alignment system
WIND    Wind Acquisition System
```

```
int Licel_TCPIP_GetCapabilities(SOCKET s, char* cap, int bufferLength);
```

```
Public Function Licel_TCPIP_GetCapabilities(ByVal client As TcpClient,
ByRef cap() As String, ByVal maxLength As Integer, ByRef validcap As Integer)
As Integer
```

Licel_TCPIP_GetMilliSecs Requests the millisecond timer value of the controller

```
int Licel_TCPIP_GetMilliSecs(SOCKET s, unsigned long* msSinceStart);
```

Transient recorder

Licel_TCPIP_SelectTR Select a Transient Recorder for subsequent communication, this command needs to send once before communicating with an individual TR like the HS_Licel_Set_Discriminator_Level and Licel_TCPIP_SetInputRange.

```
int Licel_TCPIP_SelectTR(s, int TR);
```

```
Public Function Licel_TCPIP_SelectTR(ByVal client As TcpClient,
ByVal TR As Integer) As Integer
```

Licel_TCPIP_SelectMultipleTR Select a list of Transient Recorders for subsequent communication, this needs to be send before calling the multiple routines like GetMultipleShots, GetMultipleShotsAB, Licel_TCPIP_MultipleContinueAcquisition, Licel_TCPIP_MultipleStopAcquisition.

```
int Licel_TCPIP_SelectMultipleTR(SOCKET s, int* TRList, int trNumber);
```

```
Public Function Licel_TCPIP_SelectMultipleTR(ByVal client As TcpClient,
ByVal TRList() As Integer, ByVal trNumber As Integer) As Integer
```

SetDiscriminatorLevel Set the discriminator level between 0 and 63

```
int HS_Licel_Set_Discriminator_Level(short int iDevice ,int iDiscrLevel);
```

```
Public Function HS_Licel_Set_Discriminator_Level(iDevice As Integer,
iDiscrLevel As Integer) As Integer
```

```
int Licel_TCPIP_SetDiscriminatorLevel(SOCKET s, int iDiscrLevel);
```

```
Public Function Licel_TCPIP_SetDiscriminatorLevel(ByVal client As TcpClient,
ByVal iDiscrLevel As Integer) As Integer
```


SetRange Change the input voltage range (20,100,500mV)

```
int HS_Licel_Set_Range(iDevice,int iRange);
Public Function HS_Licel_Set_Range(iDevice As Integer, iRange As Integer) As Integer
int Licel_TCPIP_SetInputRange(SOCKET s, iRange);
Public Function Licel_TCPIP_SetInputRange(ByVal client As TcpClient,
ByVal iRange As Integer) As Integer
```

SetThresholdMode Set the scale of the discriminator level. In the low threshold mode the discriminator level 63 corresponds to 25mV while in the high threshold mode it corresponds to 100mV.

```
int HS_Licel_Set_ThresholdMode(short int iDevice, int iMode);
Public Function HS_Licel_Set_ThresholdMode(iDevice As Integer, iMode As Integer) As Integer
int Licel_TCPIP_SetThresholdMode(SOCKET s, int iMode);
Public Function Licel_TCPIP_SetThresholdMode(ByVal client As TcpClient,
ByVal iMode As Integer) As Integer
```

Licel_TCPIP_GetTRTYPE Get various Transient Recorder parameter specification information

```
int Licel_TCPIP_GetTRTYPE(SOCKET s, unsigned long int * trfifoLength, int *
trSernum, int * trPCbits, int * trADCBits, double * trBinWidth
, unsigned long long* HWCAP, double* binShift);
```

Licel_TCPIP_GetFreqDivider Retrieve the frequency divider exponent, the values are valid only for units supporting this feature. To get the actual bin width multiply the bin width returned by TRTYPE with the `2FreqDivExponent`.

If the `FreqDivExponent` is 0 you have a 40MHz sampling frequency for a TR40-16bit-3U, if its is one - 20MHz, for 2 it will be 10MHz and so on.

This is the parameter that the TCPIP `FREQDIV` command accepts.

```
int Licel_TCPIP_GetFreqDivider(SOCKET s, int* freqDivExponent);
```

Licel_TCPIP_SetFreqDivider Set the frequency divider exponent, this will have effect only on units supporting this feature, it changes the sampling rate before the summation of the data.

If the `FreqDivExponent` is 0 you have a 40MHz sampling frequency for a TR40-16bit-3U, if its is one - 20MHz, for 2 it will be 10MHz and so on.

```
int Licel_TCPIP_SetFreqDivider(SOCKET s, int freqDivider);
```

Licel_TCPIP_SetShotLimit Enable or Disable the 64k Shot acquisition mode

```
int Licel_TCPIP_SetShotLimit(SOCKET s, int mode);
```

Licel_TCPIP_BlockRackTrigger Set which rack trigger should be blocked (note: the front trigger can't be blocked) Mode 0 unblocks the TR completely, this is also the default state. Mode 1 blocks trigger A, 2 blocks B... The command will only work if the hardware capabilities of the transient recorder support this. This is indicated by bit 5 in the `HWCAP` field of the `TRTYPE?` command.

```
int Licel_TCPIP_BlockRackTrigger(SOCKET s, int mode);
```

Licel_TCPIP_PreTrigger Enable the pretrigger for a selected TR. In TR20-16bit this will be 128 bins long shipped till 2018, since 2018 the TR40-16bit-3U will have 1/16 of the trace length. This means for a 16k, the pretrigger will be 1024 bins long. The TR will power up with pretrigger off. TR devices supporting pretrigger indicate it by bit 3 in the `HWCAP` field of the `TRTYPE?` command.

```
int Licel_TCPIP_PreTrigger(SOCKET s, int enablePretrigger);
```

Licel_TCPIP_SetSlaveMode Set the slave mode

```
int Licel_TCPIP_SetSlaveMode(SOCKET s);

Public Function Licel_TCPIP_SetSlaveMode(ByVal client As TcpClient)
As Integer
```

Licel_TCPIP_SetPushMode Activate the push mode for the currently selected Transient Recorder

```
int Licel_TCPIP_SetPushMode(SOCKET s, int shots,int dataType,
int numberToRead, int memory);

Public Function Licel_TCPIP_SetPushMode(ByVal client As TcpClient,
ByVal shots As Integer, ByVal dataType As Integer, ByVal numberToRead
As Integer, ByVal memory As Integer) As Integer
```

Licel_TCPIP_MPushStart Starting the MPush Mode Acquisition, this is more versatile than Licel_TCPIP_SetPushMode and should be used in new programs. See the mpush_start for a usage of this function.

```
int Licel_TCPIP_MPushStart(SOCKET s, int iShots, int* TRList, int* dataType,
int* numberToRead, int* Memories, int numPushDataSets);
```

Licel_TCPIP_GetStatus Return the status information for one Transient Recorder

```
int HS_Licel_Get_Status(short int iDevice,long int * iCycles, int * iMemory,
int * iAcq_State, int * iRecording);

Public Function HS_Licel_Get_Status(iDevice As Integer, iCycles As Integer,
iMemory As Integer, iAcq_State As Integer, iRecording As Integer) As Integer

int Licel_TCPIP_GetStatus(SOCKET s, long int* iCycles, int* iMemory,
int* iAcq_State, int* iRecording);

Public Function Licel_TCPIP_GetStatus(ByVal client As TcpClient,
ByRef shotNumber As Integer, ByRef lastMemory As Integer,
ByRef acquisitionState As Integer, ByRef recording As Integer) As Integer
```

GetMultipleShots This is the version of Licel_TCPIP_GetStatus for multiple transient recorders. It will return for each TR the number of shots accumulated over the memories and a number indicating in which memory the last shot was acquired. Newer TR supporting the separate shot counters are better served with the Licel_TCPIP_GetMultipleShotsAB routine.

```
int Licel_TCPIP_GetMultipleShots(SOCKET s, long int* shotNumbers, int maxNum,
int* numElements);
```

GetShotsAB This is returns the shotnumbers for the individual memories of the TR, if the TR does support this.

```
int Licel_TCPIP_GetShotsAB(SOCKET s, long int* shotNumbers, int maxNum,
int* numElements);
```

GetMultipleShotsAB This is the version of Licel_TCPIP_GetShotsAB for multiple transient recorders.

```
Licel_TCPIP_GetMultipleShotsAB(SOCKET s, long int* shotNumbers, int maxNum,
int* numElements);
```

ContinueAcquisition Continue the recording process without a new initialisation of the memory

```
int HS_Licel_Continue_Acquisition(short int iDevice);

Public Function HS_Licel_Continue_Acquisition(iDevice As Integer) As Integer

int Licel_TCPIP_ContinueAcquisition(SOCKET s);

Public Function Licel_TCPIP_ContinueAcquisition(ByVal client As TcpClient)
As Integer
```

Licel_TCPIP_MultipleContinueAcquisition Continue the recording process for the previously selected Transient Recorders without a new initialisation of the memory

```
int Licel_TCPIP_MultipleContinueAcquisition(SOCKET s);

Public Function Licel_TCPIP_MultipleContinueAcquisition(ByVal client
As TcpClient) As Integer
```

StopAcquisition Stop the recorder after the next received trigger.

```
int HS_Licel_Stop_Acquisition(short int iDevice);

Public Function HS_Licel_Stop_Acquisition(iDevice As Integer)As Integer

int Licel_TCPIP_Stop(SOCKET s);

Public Function Licel_TCPIP_StopAcquisition(ByVal client As TcpClient)
As Integer
```

Licel_TCPIP_MultipleStopAcquisition Stop the acquisition process for the previously selected Transient Recorders with the next received trigger pulse

```
int Licel_TCPIP_MultipleStopAcquisition(SOCKET s);

Public Function Licel_TCPIP_MultipleStopAcquisition(ByVal client
As TcpClient) As Integer
```

ClearMemory Clear both memories (A and B) of the Transient Recorder

```
int HS_Licel_Clear_Memory(short int iDevice);

Public Function HS_Licel_Clear_Memory(iDevice As Integer) As Integer

int Licel_TCPIP_ClearMemory(SOCKET s);

Public Function Licel_TCPIP_ClearMemory(ByVal client As TcpClient)
As Integer
```

Licel_TCPIP_MultipleClearMemory Clear both memories (A and B) of the previously selected Transient Recorders

```
int Licel_TCPIP_MultipleClearMemory(SOCKET s);

Public Function Licel_TCPIP_MultipleClearMemory(ByVal client As TcpClient)
As Integer
```

Licel_TCPIP_SetMaxShots Set the maximum shotnumber of the TR this can be an arbitrary number between 2 and 65335, the startup default is 4096.

```
int Licel_TCPIP_SetMaxShots(SOCKET s, int maxShots);
```

Licel_TCPIP_SetMaxBins Set the maximum number of bins of the TR unit, if the Memory DIP Switch 5 is in the ON position. This allows to adapt the TR unit better to the repetition rate and trace length requirements.

```
int Licel_TCPIP_SetMaxBins(SOCKET s, int numMaxBins);
```

StartAcquisition Start the acquisition process with a new initialization of the memory

```
int HS_Licel_Start_Acquisition(short int iDevice);

Public Function HS_Licel_Start_Acquisition(iDevice As Integer) As Integer

int Licel_TCPIP_Start(SOCKET s);

Public Function Licel_TCPIP_StartAcquisition(ByVal client As TcpClient)
As Integer
```

Licel_TCPIP_MultipleStart Start the acquisition process for the previously selected Transient Recorders with a new initialization of the memory

```
int Licel_TCPIP_MultipleStartAcquisition(SOCKET s);

Public Function Licel_TCPIP_MultipleStart (ByVal client As TcpClient)
As Integer
```

SingleShot Acquire one shot

```
int HS_Licel_Single_Shot(short int iDevice);

Public Function HS_Licel_Single_Shot (iDevice As Integer) As Integer

int Licel_TCPIP_SingleShot (SOCKET s);

Public Function Licel_TCPIP_SingleShot (ByVal client As TcpClient) As Integer
```

WaitForReady Wait for the return of the Transient Recorder from the armed state. If the waiting time is longer than the time specified by delay, then the Transient Recorder will return to the idle state with the next reading of binary data.

```
int HS_Licel_Wait_For_Ready(short int iDevice, int imDelay);

Public Function HS_Licel_Wait_For_Ready (iDevice As Integer, imDelay As
Integer) As Integer

int Licel_TCPIP_WaitForReady (SOCKET s, imDelay);

Public Function Licel_TCPIP_WaitForReady (ByVal client As TcpClient,
ByVal delay As Integer) As Integer
```

Licel_TCPIP_MultipleWaitForReady Wait until all Transient Recorders return from the armed state

```
int Licel_TCPIP_MultipleWaitForReady (SOCKET s, imDelay);

Public Function Licel_TCPIP_MultipleWaitForReady (ByVal client As TcpClient,
ByVal delay As Integer) As Integer
```

Licel_TCPIP_ReadData Read binary data into a byte array. Transient recorder data is internally 16bits wide so for every data point two bytes need to be fetched.

```
int Licel_TCPIP_ReadData (SOCKET s, int numberToRead, unsigned char* data);

Public Function Licel_TCPIP_ReadData (ByVal client As TcpClient,
ByVal numberToRead As Integer, ByRef data() As Byte) As Integer
```

GetDatasets/Read16bitwide Read binary datasets from a Transient Recorder

```
int HS_Licel_Read16bit_wide (short int iDevice, int dataType, int iNumber,
int iMemory, unsigned short * uPortData);

Public Function HS_Licel_Read16bit_wide (iDevice As Integer, dataType As
Integer, iNumber As Integer, iMemory As Integer, uPortData() As Integer)

int Licel_TCPIP_GetDatasets (SOCKET s, int iDevice, int dataType, int iNumber,
int iMemory, unsigned char* data);

Public Function Licel_TCPIP_GetDatasets (ByVal client As TcpClient,
ByVal TR As Integer, ByVal dataType As Integer, ByVal numberToRead As Integer,
ByVal memory As Integer, ByRef data() As Byte) As Integer
```

RequestDataSet Requesting the raw data sets (analog LSW, analog MSW or photon counting) from * the specified device for later read

```
int Licel_TCPIP_RequestDataSet (SOCKET s, int iDevice, int dataType,
int iNumber, int iMemory);
```

4.3 licel_data

mpush_start Activate the push mode for multiple transient recorders and multiple datasets on each TR.

```
int mpush_start(SOCKET s, int iShots, DataSet* dataSets, int nDataSetNumber,
int* numBytes);
```

push_parser transfer the binary push data into the dataSets array to be later stored in data files

```
int push_parser(unsigned char* pushBuffer, int pushOffset,
unsigned short* mem_low, unsigned short* mem_high, DataSet* dataSets,
int nDataSetNumber, int numBytes, unsigned long* timeStamp, int* dataValid);
```

4.4 licel_util

Combine_Analog_Datasets Converts the LSW and the MSW read out data values from a 12 bit ADC Transient Recorder into an integer array containing the summed up analog values. The first trash element (due to the data transmission scheme) is also removed.

```
void Licel_Combine_Analog_Datasets(unsigned short * uLSW, unsigned short *
uMSW, int iNumber, unsigned long * lAccumulated, short * iClipping);
```

```
Public Sub Licel_Combine_Analog_Datasets(iLsw() As Integer, iMsw() As
Integer, iNumber As Integer, lAccumulated() As Long, iClipping() As
Integer)
```

Combine_Analog_Datasets_16_bit Converts the LSW, the MSW and the PHM read out data values from a 16 bit ADC Transient Recorder into an integer array containing the summed up analog values. The first trash element (due to the data transmission scheme) is also removed.

```
void Licel_Combine_Analog_Datasets_16_bit(unsigned short* uLSW, unsigned short*
uMSW, unsigned short* uPHM, int iNumber, unsigned long* lAccumulated,
short* iClipping);
```

Licel_Combine_Analog_Squared_Data Converts the LSW, MSW and HSW values of the squared data into an 64 bit integer array containing the summed up squared analog values. The first trash element (due to the data transmission scheme) is also removed

```
void Licel_Combine_Analog_Squared_Data(unsigned short* uSQLSW, unsigned short*
uSQMSW, unsigned short* uSQHSW, int iNumber, unsigned long long*
llSQAccumulated);
```

Convert_Photoncounting Converts 16 bits of raw Photon counting data into an integer array containing the summed up photon counting values. The first trash element (due to the data transmission scheme) is also removed. The clipping information present in the most significant bit is masked out if necessary*/

```
void Licel_Convert_Photoncounting(unsigned short * photon_raw, int iNumber,
unsigned long * photon_c, int iPurePhoton);
```

```
Public Sub Licel_Convert_Photoncounting(photon_raw() As Integer, iNumber As
Integer, photon_c() As Long, iPurePhoton As Boolean)
```

Convert_Photoncounting_FullWord Converts 24 bits of raw Photon counting data into an integer array containing the summed up photon counting values. The first trash element (due to the data transmission scheme) is also removed.

```
void Licel_Convert_Photoncounting_FullWord(unsigned short* uPHO,
unsigned short* uPHM, int iNumber, unsigned long* photon_c);
```

Licel_Combine_Photon_Squared_Data Converts the LSW, MSW values into an 64 bit integer array containing the summed up squared photon counting values. The first trash element (due to the data transmission scheme) is also removed.

```
void Licel_Combine_Photon_Squared_Data(unsigned short* uSQLSW,
unsigned short* uSQMSW, int iNumber, unsigned long long* llSQAccumulated);
```

Licel_GetSquareRootBinary Convert the squared data to binary number for the std dev

```
void Licel_GetSquareRootBinary(unsigned long* lAccumulated,
unsigned long long* llSQAccumulated, int iNumber, int iShots,
unsigned long* sqd_bin);
```

Normalize_Data Normalizes the accumulated Data with respect to the number of cycles

```
void Licel_Normalize_Data( unsigned long * lAccumulated, int iNumber,
int iShots, double * dNormalized);
```

```
Public Sub Licel_Normalize_Data(lAccumulated() As Long, iNumber As Integer,
iShots As Integer, dNormalized() As Double)
```

Licel_Normalize_SquaredData Normalizes the accumulated squared Data with respect to the number of shots

```
void Licel_Normalize_SquaredData(unsigned long* sqd_bin, int iNumber,
int iShots, double* dSampleStandardDev);
```

Licel_MeanError convert the sample standard deviation to the more relevant error of the mean value

```
void Licel_MeanError(double* dSampleStandardDev, int iNumber, int iShots,
double* meanError);
```

Scale_Analog_Data Scales the normalized data with respect to the input range

```
void Licel_Scale_Analog_Data(double * dNormalized, int iNumber, int
iRange, int trADCBits, double * dmVData);
```

```
Public Sub Licel_Scale_Analog_Data(dNormalized() As Double, iNumber As
Integer, iRange As Integer)
```

ERRMess Error message dump with error code and where the error occurred

```
void ERRMess(int nNumber, const char *sPlace);
```

4.5 APD functions

Licel_TCPIP_APDGetStatus Get the status of the APD with the corresponding APD number.

```
int Licel_TCPIP_APDGetStatus(SOCKET s, int APD, bool* ThermoCooler,
bool* TempInRange, int* HV, bool* HVControl);
```

```
Public Function Licel_TCPIP_APDGetStatus(ByVal client As TcpClient,
ByVal APD As Integer, ByVal ThermoCooler As Boolean, ByRef TempInRange
As Boolean, ByRef HV As Integer, ByRef HVControl As Boolean) As Integer
```

Licel_TCPIP_APDSetCoolingState Set the cooling mode of the specified APD

```
int Licel_TCPIP_APDSetCoolingState(SOCKET s, int APD, bool ThermoCooler);
```

```
Public Function Licel_TCPIP_APDSetCoolingState(ByVal client As TcpClient,
ByVal APD As Integer, ByVal ThermoCooler As Boolean) As Integer
```

Licel_TCPIP_APDSetGain Set the applied high voltage (gain) of the specified APD

```
int Licel_TCPIP_APDSetGain(SOCKET s, int APD, int HV);
```

```
Public Function Licel_TCPIP_APDSetGain(ByVal client As TcpClient,
ByVal APD As Integer, ByVal HV As Integer) As Integer
```

4.6 PMT functions

Licel_TCPIP_PMTGetStatus Get the status of the PMT with the corresponding PMT number

```
int Licel_TCPIP_PMTGetStatus(SOCKET s, int PMT, bool* HVOn, float* HV,
bool* HVControl);
```

```
Public Function Licel_TCPIP_PMTGetStatus(ByVal client As TcpClient,
ByVal PMT As Integer, ByRef HVOn As Boolean, ByRef HV As Double,
ByRef HVControl As Boolean) As Integer
```

Licel_TCPIP_PMTSetGain Set the applied high voltage (gain) of the specified PMT

```
int Licel_TCPIP_PMTSetGain(SOCKET s, int PMT, int HV);
```

```
Public Function Licel_TCPIP_PMTSetGain(ByVal client As TcpClient,
ByVal PMT As Integer, ByVal HV As Integer) As Integer
```

4.7 Timing functions

The old API was designed for a single trigger board, if one Ethernet Controller controls more than one trigger board, each board needs to be addressed separately with a board ID. The first board has the default ID 0, which is addressed by the old API, the additional boards have the ID's 1 and 2.

Licel_TCPIP_SetTriggerMode Enable/Disable the trigger in and outputs.

Old API defaults to **boardID**: 0.

```
int Licel_TCPIP_SetTriggerMode(SOCKET s, bool LaserActive,
bool PreTriggerActive, bool QSwitchActive, bool GatingActive,
bool MasterTrigger);
```

```
Public Function Licel_TCPIP_SetTriggerMode(ByVal client As TcpClient,
ByVal LaserActive As Boolean, ByVal PretriggerActive As Boolean,
ByVal QSwitchActive As Boolean, ByVal GatingActive As Boolean,
ByVal MasterTrigger As Boolean) As Integer
```

New API for multiple trigger boards.

```
int Licel_TCPIP_SetTriggerMode(SOCKET s, int boardID, bool LaserActive,
bool PreTriggerActive, bool QSwitchActive, bool GatingActive,
bool MasterTrigger);
```

```
Public Function Licel_TCPIP_SetTriggerModeN(ByVal client As TcpClient,
ByVal boardID As Integer, ByVal LaserActive As Boolean,
ByVal PretriggerActive As Boolean, ByVal QSwitchActive As Boolean,
ByVal GatingActive As Boolean, ByVal MasterTrigger As Boolean) As Integer
```

Licel_TCPIP_SetTriggerTiming Set the timing parameter, as for the trigger mode there is the old API, which defaults to **boardID**: 0 while the new API supports multiple trigger boards.

Old API:

```
int Licel_TCPIP_SetTriggerTiming(SOCKET s, long repetitionRate,
long Pretrigger, long PretriggerLength, long QSwitch, long QswitchLength);
```

```
Public Function Licel_TCPIP_SetTriggerTiming(ByVal client As TcpClient,
ByVal repetitionRate As Integer, ByVal Pretrigger As Integer,
ByVal PretriggerLength As Integer, ByVal QSwitch As Integer,
ByVal QswitchLength As Integer)
```

New API:

```
int Licel_TCPIP_SetTriggerTiming(SOCKET s, int boardID, long repetitionRate,
long Pretrigger, long PretriggerLength, long QSwitch, long QswitchLength);
```



```
Public Function Licel_TCPIP_SetTriggerTiming(ByVal client As TcpClient,
ByVal boardID As Integer, ByVal repetitionRate As Integer, ByVal Pretrigger
As Integer, ByVal PretriggerLength As Integer, ByVal QSwitch As Integer,
ByVal QswitchLength As Integer)
```

4.8 Security functions

Licel_TCPIP_SetAccessLimited Activate the access limitation, that means only whitelisted hosts can access the Controller and need to verify them self by properly encoding with the connectionPasswd a two 8 byte numbers. Make sure that you called [Licel_TCPIP_SetWhiteList](#) before, otherwise no host will be authorized to access the Controller.

```
int Licel_TCPIP_SetAccessLimited(SOCKET s, char* connectionPasswd, char* passwd);
```

Licel_TCPIP_SetAccessUnLimited Deactivate the access limitation, that means every hosts can access the Controller.

```
int Licel_TCPIP_SetAccessUnLimited(SOCKET s, char* passwd);
```

Licel_TCPIP_SetWhiteList List hosts that be allowed to to access the Controller after [Licel_TCPIP_SetAccessLimited](#) has been called. One can list three different hosts. Specifying a 255 as the last number activates the whole range, e.g. 10.49.234.255 as host will make all hosts from 10.49.234.1 to 10.49.234.254 whitelisted hosts.

```
int Licel_TCPIP_SetWhiteList(SOCKET s, char* whiteHost1, char* whiteHost2,
char* whiteHost3, char* passwd);
```

4.9 Power Meter

The power meter uses two data sockets the first socket for command transmission and return values. The second which has port number one above the first for continuous data transmission once the data acquisition has been started.

Licel_TCPIP_PowerSelectChannel Select one channel of the power meter.

```
int Licel_TCPIP_PowerSelectChannel(SOCKET s, int Channel );
```

Licel_TCPIP_PowerStart Start the power meter data acquisition.

```
int Licel_TCPIP_PowerStart(SOCKET s);
```

Licel_TCPIP_PowerTrace Start a single power meter data acquisition, it will return a raw integer array of ADC readings.

```
int Licel_TCPIP_PowerTrace(SOCKET s, int *readings, int *numReadings);
```

Licel_TCPIP_PowerStop Stop the power meter data acquisition.

```
int Licel_TCPIP_PowerStop(SOCKET s);
```

Licel_TCPIP_PowerGetData Get the data from the second socket with the port number which is one above the command socket.

```
int Licel_TCPIP_PowerGetData(SOCKET s, int *milliSeconds, double *reading);
```

4.10 Bore Alignment

The bore alignment detector uses two data sockets the first socket for command transmission and return values. The second which has port number one above the first for continuous data transmission once the data acquisition has been started.

Licel_TCPIP_BoreSetRanges Set the background and the signal region for the bore alignment sensor.

```
int Licel_TCPIP_BoreSetRanges(SOCKET s, int backgroundStart, int backgroundStop,
int signalStart, int signalStop);
```

Licel_TCPIP_BoreSign Toggle the sign of the counter that is transmitted together with the alignment data. The counter will increment with every cycle, however the sign might toggle. This can be used to make sure that the data evaluated has been measured after a certain point, for instance a drive movement.

```
int Licel_TCPIP_BoreSign(SOCKET s);
```

Licel_TCPIP_BoreStart Start the alignment sensor data acquisition. The data itself will be transmitted over a second socket (see [BoreGetData](#))

```
int Licel_TCPIP_BoreStart(SOCKET s, int shots, int cycles);
```

Licel_TCPIP_BoreStop Stop the alignment sensor data acquisition.

```
int Licel_TCPIP_BoreStop(SOCKET s);
```

Licel_TCPIP_BoreGetData Get the data from the second socket with the port number which is one above the command socket. The data is returned in the [countrates](#) array which should hold 8 doubles.

```
int Licel_TCPIP_BoreGetData(SOCKET s, double *countrates, long int *counter);
```

4.11 Wind Acquisition

The wind acquisition system uses two data sockets the first socket for command transmission and return values. The second which has port number one above the first for continuous data transmission once the data acquisition has been started.

Licel_TCPIP_WindSetShots Set the shot number

```
int Licel_TCPIP_WindSetShots(SOCKET s, int shots);
```

Licel_TCPIP_WindSetFFTSize Set the number of ADC Samples per FFT

```
int Licel_TCPIP_WindSetFFTSize(SOCKET s, int FFTsize);
```

Licel_TCPIP_WindSetRuns Set the number of runs to be acquired

```
int Licel_TCPIP_WindSetRuns(SOCKET s, int runs);
```

Licel_TCPIP_WindSetNumFFT Set the number of runs to be acquired

```
int Licel_TCPIP_WindSetNumFFT(SOCKET s, int numFFT);
```

Licel_TCPIP_WindSetParameter Send WIND acquisition parameters via TCPIP.

```
int Licel_TCPIP_WindSetParameter(SOCKET s, int shots, int FFTsize, int runs,
int numFFT);
```

Licel_WindCalcRangeBins Calculate the number of FFTs resulting by choosing the distance and FFT Size.

```
int Licel_WindCalcRangeBins(int distance, int FFTsize, int *numBytes);
```

Licel_TCPIP_WindAcqStart Start the acquisition.

```
int Licel_TCPIP_WindAcqStart(SOCKET s);
```

Licel_WindPushParser Transfer the binary push data into the [dataSets](#) array to be later stored in data files.

```
int Licel_WindPushParser(unsigned char* pushBuffer, int pushOffset, int numBytes,
int* dataValid, unsigned long long* Data, unsigned long* timeStamp);
```

4.12 Function arguments

<i>trfifoLength</i>	<i>Indicates the set FIFO Length of the Transient Recorder.</i>
<i>trSerNum</i>	<i>Indicates the serial number of the Transient Recorder.</i>
<i>trPCbits</i>	<i>Indicates the number of Photon Counting Bits.</i>
<i>trADCBits</i>	<i>Indicates the number of ADC bits of the Transient Recorder. This can be 12 or 16 bits.</i>
<i>trBinWidth</i>	<i>Indicates the base width in m of a single used by the Transient Recorder for ADC and the photon counting</i>
<i>HWCAP</i>	<p><i>Further hardware capabilities of the transient recorder, this is a bitfield with the following bits currently defined.</i></p> <p><i>0x01 separate shot counter B</i></p> <p><i>0x02 separate shot counter C</i></p> <p><i>0x04 separate shot counter D</i></p> <p><i>0x08 pretrigger</i></p> <p><i>0x10 memory blocking</i></p> <p><i>0x20 squared data support</i></p> <p><i>0x40 freq divider</i></p>
<i>binShift</i>	<i>bin shift of the analog data with respect to the PC data</i>
<i>freqDivider</i>	<i>if the HWCAP bit 0x40 is set the TR supports changing the frequency. The bin width reported by TRType needs then to be multiplied with freqDivider. Possible values are 1, 2, 4, 8, 16, 32, 64, 128.</i>
<i>mode</i>	<i>Mode 0 - turn off the 64k shot acquisition capability. Mode 1 - turn on the 64k shot acquisition capability.</i>
<i>enablePretrigger</i>	<i>enablePretrigger 0 - turn off the pretrigger trace section, enablePretrigger 1 - turn on the pretrigger trace.</i>
<i>shotNumber</i>	<i>Number of shots already acquired. This shot number has an offset of 2 as the two initial clearing cycles advance the shot number to 2</i>
<i>shotNumbers</i>	<i>Array with a number of shots already acquired. If the AB variants are used it will list the number of shots per memory. For MSHOTS it will show the cumulative number of shots for both memories as shotnumber above and a 0 if the previous shot was on memory A and 1 if it was on memory B.</i>
<i>maxNum</i>	<i>Maximum number of elements the shotnumbers array can hold.</i>
<i>numElements</i>	<i>Number of elements set in shotnumbers by the subroutine.</i>
<i>maxShots</i>	<i>Maximum number of shots before the TR stops itself without a stop command. The default at start up is 4096.</i>
<i>numMaxBins</i>	<i>define the maximum number of bins in the memory. Valid range 2 to 32768. The numberToRead in Licel.TCPIP_GetDataset should be smaller than this number otherwise invalid data will be returned for all bins where the index is larger than numMaxBins.</i>
<i>lastMemory</i>	<i>Memory to which the previous acquisition was added.</i>
<i>acquisitionState</i>	<i>FALSE when the transient returns from the armed state, TRUE when an acquisition is running</i>
<i>recording</i>	<i>TRUE during acquisition-time, e.g. the ADC or the photon counting is acquiring data. FALSE during summation and when the TR is waiting for a new trigger.</i>
<i>numberToRead</i>	<i>Number of 16 bit wide data points</i>

<i>dataType</i>	States which part of the raw information should be transferred from the device to the computer. Use the constants PHOTON, LSW, MSW and PHM !!! Note: PRxx xx recorders need to read LSW and MSW instead of PHOTON, for squared analog data A2L, A2M and A2M should be used. And for squared photon counting data P2L and P2M shall be used. Use the following symbolic constants PHOTON 0 LSW 1 MSW 2 PHM 3 P2L 4 P2M 5 A2L 6 A2M 7 A2H 8
<i>dataSets</i>	Array of structures which hold the information about the datasets like number of bins, Range, discriminator threshold etc.
<i>nDataSetNumber</i>	number of elements in the dataSets array.
<i>numPushDataSets</i>	number of push datasets for the MPUSH Mode
<i>Memories</i>	list of summation memories to be retrieved in the mpush mode.
<i>numBytes</i>	number of bytes in the mpush data transmission blocks.
<i>pushBuffer</i>	the push buffer containing the raw bytes
<i>pushOffset</i>	how many bytes are in the pushBuffer
<i>mem_high</i>	buffer for the MSW data (U16)
<i>mem_low</i>	buffer for the LSW data (U16)
<i>dataValid</i>	1 if the mpush data data has been transferred into the dataSets array, 0 otherwise.
<i>timeStamp</i>	time stamp for the mpush data, it will only be updated if dataValid is 1. The value is the number of milliseconds since controller start
<i>memory</i>	Summation memory to be retrieved.
<i>shots</i>	Transfer data every #shots.
<i>TR</i>	Hardware addresses of the Transient Recorder. Valid values are 0:15. All single device commands to the Ethernet Controller will access the corresponding Transient Recorder selected.
<i>TRLList</i>	List containing the hardware addresses of the Transient Recorders for subsequent "multiple" commands. All multiple device commands to the Ethernet Controller will access the Transient Recorders mentioned in this list.
<i>trNumber</i>	Length of the TRLList.
<i>iDevice</i>	Hardware address of the Transient Recorder when used with the DIO-32HS card. Valid values are 0:7.
<i>s, client</i>	TCPIP socket connection reference
<i>iDiscrLevel</i>	Photon counting discriminator level. Valid values are 0:63
<i>iRange</i>	Analog input range valid values 0:2. Please use symbolic constants, as defined below. MILLIVOLT500 0 MILLIVOLT100 1 MILLIVOLT20 2


<i>iMode</i>	Photon counting threshold mode valid values 0:1. Please use symbolic constants, as defined below. THRESHOLD_LOW 1 THRESHOLD_HIGH 0
<i>iCycles</i>	Number of already acquired traces. Please note that there are two cycles used for memory initialization, so that the current shot number is <i>iCycle</i> -2 if <i>iCycle</i> >=2. Valid values are: 0:4095
<i>iShots</i>	Number of acquired shots. This is <i>iCycles</i> - number of memories in the TR
<i>iMemory</i>	Indicates to which memory bank the last shot was added. Valid values are 0:1
<i>iAcq_State</i>	0 if there is an acquisition currently running, otherwise 1 .
<i>iRecording</i>	1 if the ADC is acquiring values, when the status command was issued, otherwise 0.
<i>imDelay</i>	Delay to wait (in milliseconds).
<i>iNumber</i>	Number of Data points. The maximum value depends on the type of the Transient Recorder. for TR xx:80 8192 for TR xx:160 16380
<i>uPortData</i>	Array to store the data values.
<i>uLSW</i>	Array containing the LSW readout value of the analog data.
<i>uPHO</i>	Array containing the PHO readout value of the photon counting data.
<i>uMSW</i>	Array containing the MSW readout value of the analog data.
<i>uPHM</i>	Array containing the PHM readout value of the photon counting/analog data.
<i>uSQLSW</i>	Array containing the A2L (LSW) readout value of the squared analog data.
<i>uSQMSW</i>	Array containing the A2M (MSW) readout value of the squared analog data.
<i>uSQHSW</i>	Array containing the A2H (HSW) readout value of the squared analog data
<i>IAccumulated</i>	Array containing the summed up analog data.
<i>IISQAccumulated</i>	Array containing the summed up squared analog data.
<i>sqd_bin</i>	square root of (<i>iShots</i> * <i>IISQAccumulated</i> - <i>IAccumulated</i> ²) see eq. 11.
<i>dSampleStandardDev</i>	sample standard deviation for each bin see eq. 12
<i>meanError</i>	mean value standard deviation for each bin see eq. 13
<i>iClipping</i>	Array containing the clipping(out of range) information. 1 if the overange condition (for the specific data point) is at least fulfilled once, otherwise 0.
<i>photon_raw</i>	Array containing the raw photon counting data.
<i>photon_c</i>	Array with photon counting data without clipping information.
<i>iPurePhoton</i>	TRUE if there is no analog data (hence no need to remove the clipping bit). FALSE otherwise.
<i>dNormalized</i>	Contains the array normalized to the shot number.
<i>dmVData</i>	Array converted to mV.
<i>sHost</i>	String with the host name.
<i>iPort</i>	Integer port number to connect with the host.
<i>silent</i>	Boolean suppress MsgBox when the connection fails.

<i>command</i>	<i>String that will be transferred to the Ethernet Controller.</i>
<i>response</i>	<i>String containing the response from the Controller, the trailing CRLF will be removed.</i>
<i>maxlength</i>	<i>Max storage capacity of the response string.</i>
<i>nTimeoutMillisec</i>	<i>Max time to wait for a closing CRLF.</i>
<i>Data</i>	<i>Byte array to store the data.</i>
<i>points</i>	<i>Number of bytes to read.</i>
<i>nTimeOutMillisec</i>	<i>Max time to wait to read the specified amount of data.</i>
<i>newPort</i>	<i>New Port to connect too after reboot of the Controller.</i>
<i>passwd</i>	<i>String containing the current password for the Controller.</i>
<i>newHost</i>	<i>String containing the IP address that the Controller should be set to.</i>
<i>mask</i>	<i>String with the subnet mask that the Controller should use for TCP/IP communication.</i>
<i>gateway</i>	<i>String with the gateway that should be used by the Controller for TCP/IP communication.</i>
<i>buffer</i>	<i>String to hold the identification information.</i>
<i>bufferLength</i>	<i>Max capacity of the result string.</i>
<i>cap</i>	<i>Array containing the information about the available capabilities.</i>
<i>maxLength</i>	<i>Max. capacity of the array.</i>
<i>validcap</i>	<i>Number of different capabilities.</i>
<i>delay</i>	<i>Max. time to wait to return to the idle state (in milli seconds).</i>
<i>msSinceStart</i>	<i>milliseconds since start of the controller</i>
<i>APD</i>	<i>The physical device number of the APD. Valid values are 0:3.</i>
<i>TempInRange</i>	<i>TRUE if the thermocooler is on and the detector temperature is very close to the target temperature.</i>
<i>HVControl</i>	<i>TRUE if remote control of the high voltage is active.</i>
<i>HV</i>	<i>The applied HV to the detector.</i>
<i>ThermoCooler</i>	<i>Turns the thermo cooler ON, if FALSE the detector is only passively cooled.</i>
<i>PMT</i>	<i>The physical device number of the PMT. Valid values are 0:7.</i>
<i>HVOn</i>	<i>TRUE is the high voltage is ON.</i>
<i>remote</i>	<i>TRUE if remote control is active.</i>
<i>boardID</i>	<i>0 for the first timing board (default), 1 and 2 for the additional timing boards.</i>
<i>LaserActive</i>	<i>If TRUE a trigger for the laser lamp will be generated.</i>
<i>PretriggerActive</i>	<i>If TRUE a trigger for the Transient Recorder will be generated.</i>
<i>QSwitchActive</i>	<i>If TRUE a trigger for the laser Q-Switch will be generated.</i>
<i>GatingActive</i>	<i>If TRUE a gating pulse will be generated. The gating pulse starts with the raising edge of the pretrigger and ends with the falling edge of the Q-Switch Pulse.</i>
<i>MasterTrigger</i>	<i>If TRUE an external trigger will be accepted, else if FALSE the internal trigger will be used. The internal trigger will be controlled via the repetitionRate parameter.</i>
<i>repetitionRate</i>	<i>The internal mode delay between two pulses (in nano seconds).</i>


<i>Pretrigger</i>	<i>Delay between internal or external trigger and pre-trigger (in nano seconds).</i>
<i>PretriggerLength</i>	<i>Length of the pre-trigger pulse (in nano seconds).</i>
<i>QSwitch</i>	<i>Delay between pre-trigger start and Q-Switch start (in nano seconds).</i>
<i>QswitchLength</i>	<i>Length of the Q-Switch pulse (in nano seconds).</i>
<i>whiteHost</i>	<i>Host that is allowed to open a connection to the Controller in the limited access mode. Specifying a 255 as the last number activates a IP range, e.g. 10.49.234.255, as host will make all hosts from 10.49.234.1 to 10.49.234.254 whitelisted hosts.</i>
<i>connectionPasswd</i>	<i>Password used for encrypting the tokens sent from the Controller initially.</i>
<i>Channel</i>	<i>Power meter detector channel, valid values are 0:3.</i>
<i>readings</i>	<i>Power meter raw data for a single trace.</i>
<i>numReadings</i>	<i>Number of valid data points in the power raw trace.</i>
<i>milliSeconds</i>	<i>Milliseconds since start.</i>
<i>reading</i>	<i>Power meter reading.</i>
<i>backgroundStart</i>	<i>First background bin.</i>
<i>backgroundStop</i>	<i>Last background bin.</i>
<i>signalStart</i>	<i>First signal bin.</i>
<i>signalStop</i>	<i>Last signal bin.</i>
<i>cycles</i>	<i>Number of cycles (data transmissions), -1 for infinite cycles.</i>
<i>countrates</i>	<i>Array with 8 count rates (4 signal and 4 background).</i>
<i>counter</i>	<i>Number of transmitted packages, the sign might be negative.</i>
<i>FFTsize</i>	<i>number of points per spectra.</i>
<i>runs</i>	<i>number of runs, each run contains shots acquisitions.</i>
<i>numFFT</i>	<i>number of FFTs resulting by choosing the distance and FFT Size.</i>
<i>distance</i>	<i>max. trace distance in meter.</i>

4.13 Timing Parameter explanation

External trigger MasterTrigger = True



Internal trigger MasterTrigger = False



The Laser Lamp pulse has a fixed length of $5\mu\text{s}$.

5 Transient recorder

5.1 Low Level Commands

This applies to the parallel bus communication with the DIO32 cards and the communication of the Ethernet Controller itself with the Transient Recorder. The commands here indicate the low level operating mechanism of the Transient Recorder.

The commands are transmitted over port 2 and 3 (group2) without double buffering, while the information from the Transient Recorder is transmitted over port 0 and 1 (group 1) with double buffering. Each command also transfers a 16bit value from the device to the computer. The handshake lines of group 1 and 2 are connected. In this way a reading operation also sends a command to the device. For instance, in order to get the status information, first the corresponding command is written to port where the output signal levels remain unchanged until the next command is sent. Thus, for a sequence of reads, the output levels do not change. The commands consist of two bytes in which the first three bits of the first byte (Port 2) and the first bit of the second byte (Port 3) are the address. Bits 3-7 of the first byte contain the task information. These tasks are described below.

<i>0: Reset</i>	<i>The shot number is set to 0, the memory is not touched, the device is in the idle state.</i>
<i>8: Start</i>	<i>Used in to clear the memory and to restart an acquisition. Before a new acquisition can be started after a Reset, the two summation memories should be cleared. This is done by the Sequence 0, 24, 8, then the memory A is cleared. After the clearing a sequence 24+128,8 clears Memory B and the device is in the armed state. As an alternative, the sequence 24+128, 8, 16 returns the device into the idle state and can be restarted with an 8. An external trigger is not necessary for clearing the memory since an internal one is used. Thus clearing can be understood as an internal 0 shot.</i>
<i>16: Stop</i>	<i>Stops the acquisition process with the next received trigger, the device returns to the idle state. The shot number is not influenced by this command.</i>
<i>24: Set Memory</i>	<i>Selects the memory to be addressed (A), 24+128 - memory B.</i>
<i>32: Select DataSet</i>	<i>After the 104 command the selected data set is Photon counting. 32 selects the Analog LSW, and 32+128 selects the Analog MSW.</i>
<i>40: Advance Bin</i>	<i>The next binary value of the specified data set is sent to the PC.</i>
<i>48: Set Range</i>	<i>The Input Range is set to -100mV. Using 48+128 the Input Range is set to -20mV.</i>
<i>56 : Reset Range</i>	<i>The Input Range is set to -500mV.</i>
<i>64: Decrease Discriminator Level</i>	<i>Decreases the Discriminator Level by one.</i>
<i>72: ResetDiscriminator Level</i>	<i>The Discriminator level is reset to 63.</i>
<i>80: Status</i>	<i>Returns the status information. The status information holds the shot number in bit 0-12. The information as to whether there is an actual recording process is stored in bit 13(1=active). The information in to which memory the last shot was added is stored in 14 (0=Memory B, 1=Memory A). In bit 15 the acquisition state is stored (0=idle, 1=armed).</i>
<i>88: Reset Damping</i>	<i>The device is in the low threshold mode.</i>
<i>96: High Threshold Mode</i>	<i>The device is in the high threshold mode.</i>
<i>104: Reset Acquisition State</i>	<i>The device returns to the idle state without resetting the shot number.</i>
<i>112: Set Max Shots</i>	<i>Transmit the desired shotnumber by serializing a 16 bit bitfield with the highest bit sended first. For Zero send 112 and for 1 send 112 + 128 = 240. To set the desired shotnumber to 20 send the sequence 0000000000010100.</i>

128: Group II

Used to shift the meaning of command. See for instance 48 Set Range.

208: Status + Group II

Get the Transient Recorder hardware configuration. This will work only with Transient Recorders shipped after October 2009. To distinguish between both issue a *Power reset* before. For the old TR bit 15 will be zero and for new ones, which support the command, it will be 1.

bit	meaning
15	always 1
14,13	ADC Type 00 - 12 bit, 10 - 16 bit
12	reserved
11-8	FIFO Length see the TR Manual, the DIP switches are mirrored here
7	64k can be enabled
6,5	Width of the photon counting 00 - 4bit, 01 - 6 bit, 10 -8 bit.
4	reserved
3-0	Shotnumber bits 16 -13

As an example, to stop device 2, one has to send 18=16+2 over port 2 and 0 over port 3.


5.2 Memory organization

Current Transient Recorder The Transient Recorder has two separate memories for corresponding to trigger A and B. In each of these memory regions the data is a 57 bit wide vector of accumulated values. The length of these vectors (trace) is defined by the FIFO length as set by the dip switches. The accumulated analog data for each bin can be up to 32 bits wide plus 1 extra bit, which is used for the clipping information. The photon counting takes the remaining 24 bits.


The photon counting data can be transferred with one read operation for each bin if the shot number is not larger than 4094 shots. The analog data is transferred in typically in two read operations. Only for 16 bit ADC TR's when shot numbers larger than 32767 are used, a third read access is required. The analog data has a additional flag, indicating that the sum incorporates a overflow value. If for one bin a 12 bit ADC TR gives either 0 or 0xFFF the flags is set for the sum and indicates that either an over- or underflow has occurred at this special bin. For a 16 bit ADC TR those values would be 0 and 0xFFFF. The sum at these points may not correspond to the physical mean value as the actual ADC value could have been -10 or above 4095 (65335). This flag persists when the next trace is added to the previous traces. After accumulating, for instance 4094 shots, one is able to verify that all mean values do not incorporate out of range values by checking this clip flag. The clip flag is cleared by clearing the memory. The clip flag is transmitted as the 24th bit in the analog dataset for a 12 bit ADC TR and as the 28th bit for a 16 bit ADC TR.

By default the Transient Recorder will stop acquisition after 4094 shots. The clipping bit is then just one bit above the averaged analog data. This makes the units to behave like older Transient Recorders described below. For longer acquisitions the 64k shot mode must be activated (supported only in the newer Transient Recorders). The data bits will be above clipping bit. For the 16 bit Transient Recorder in 64k shot mode the accumulation result will be 32 bit wide and the bit 31 is mapped into the photon counting most significant word as shown in the memory organization structure below.


Analog data structure - 12 bit



Analog data structure - 16 bit





Photon counting data structure




Squared Data The squared data has 48 bits for the analog squared data and 32 bits for the squared photon data. In order to work correctly with this data its mandatory to use a 64 bit unsigned int at least for the analog squared data. **Note:** To get this working correctly the transient recorder needs to have the corresponding [hardware capability](#) and the Ethernet controller software should be newer than 2019-12-17. The software version of the Ethernet controller can be retrieved with the `*IDN?` command or with the `Search controller.vi`

Analog squared data structure - 48 bit




Photon counting squared data structure - 32 bit




Pretrigger The Pretrigger mode the data from the ADC and the photon counts are buffered. The buffer is 1/16th of the trace length. If oversampling is activated it shortens the buffer by the oversampling factor.

For a 16k trace this results in a 1k pretrigger. The buffered data covers then 25ns * 1024. If oversampling by factor of 4 (10MHz) is activated the pretrigger is 100ns * 256 with the identical time as before. For TR with 40MHz, the range bin 0 will then contain the data from 25ns * 1024 before the trigger and bin 1024 will correspond to the range bin 0 without pretrigger.


without pretrigger - 40MHz



with pretrigger - 40MHz




with pretrigger - 10MHz




Transient recorder before Oct. 2009 The Transient Recorder has two separate memories corresponding to trigger A and B. In each of these memory regions the data is a 40 bit wide vector of accumulated values. The length of these vectors (trace) is defined by the FIFO length and is 8K for the TR-xx-80 and 16k for the TR-xx-160.

The accumulated analog data for each bin is 24 bits wide and the photon counting takes the remaining 16 bits. The photon counting data can be transferred with one read operation for each bin while the analog data is transferred in two read operations. The analog data has a additional flag indicating that the sum incorporates a overflow value. If for one bin the ADC gives either 0 or 0xFFFF, the flags is set for the sum and indicates that either an over- or underflow has occurred at this special bin. The sum at these points may not correspond to the physical mean value as the actual ADC value could have been -10 or above 4095. This flag persists when the next trace is added to the previous traces. After accumulating for instance 4094 shots one is able to verify that all mean values do not incorporate out of range values by checking this clip flag. The clip flag is cleared by clearing the memory. The clip flag is transmitted as the 24th bit in the analog dataset or the 15 bit in the photon counting dataset.


Analog data structure



Photon counting data structure



PRxx-xx recorders differ in the memory layout. Here the maximum number of counts per bin for a single shot is 63 which corresponds to 6 bits. Together with 4094 shots the accumulated data can be 18 bits wide. The data is then transferred in a LSW and a MSW dataset (exactly as in the case with analog data). There is no clip bit as it is useful only for analog acquisitions.



5.3 Raw Data to Physical Value Conversion

The [Licel data file format](#) stores the data as raw values and defers the computation of physical values to the display phase.

The conversion starts with a normalization with the shot number. After this step the analog data shows the mean ADC bit values, while the photon counting shows the mean counts per bin per shot (this is the data display used by the Track and Live Display VI's).

The analog data needs then to be scaled by the ADC max value and the input range.

$$phys = norm * \frac{analogRange}{2^{ADCbits} - 1} \tag{1}$$

for a 12 bit ADC and 500mV range this means

$$mVData = norm * \frac{500mV}{4095} \tag{2}$$

The photon counting data can be converted from the counts per bin per shot into MHz if number of bins per μs is given as:

$$MHzData = norm * \frac{bins}{\mu s} \tag{3}$$

If for instance the counts per bin per shots are 1.5 and the number of bins per μs is 20, this would correspond to 30MHz. The Transient Recorder units share the clock between the ADC and photon counting so the number of bins per microsecond and the sampling rate are equal.

5.4 Raw Squared data to signal standard deviation

The sum of counts c at each bin after acquiring N shots is read from the transient recorders (standard photon counting or analog memories):

$$c_{bin} = \sum_{i=0}^{N-1} x_i \quad (4)$$

with the counts x_i of the i th shot.

The mean value μ (counts per shots) acquired in a bin is then

$$\mu = \frac{c_{bin}}{N} = \frac{\sum_{i=0}^{N-1} x_i}{N}. \quad (5)$$

The *sample standard deviation* is (https://en.wikipedia.org/wiki/Standard_deviation):

$$s = \sqrt{\frac{\sum_{i=0}^{N-1} (x_i - \mu)^2}{N - 1}}. \quad (6)$$

From there we get the *standard error of the mean* (https://en.wikipedia.org/wiki/Standard_error):

$$\sigma_{\mu} \approx \frac{s}{\sqrt{N}}. \quad (7)$$

As the x_i are the counts for each shot, but the counts are already summed by the transient recorders there is no direct access to standard deviation nor to standard error when reading from standard analog or PC memories. The solution is to read squared data c_{bin}^2 from transient recorders supporting the summation of squared counts. Then one can calculate the sample standard deviation using

$$\begin{aligned} s &= \sqrt{\frac{c_{bin}^2 - (c_{bin})^2/N}{N - 1}} \\ &= \sqrt{\frac{\sum_{i=0}^{N-1} x_i^2 - (\sum_{i=0}^{N-1} x_i)^2/N}{N - 1}} \\ &= \sqrt{\frac{N \sum_{i=0}^{N-1} x_i^2 - (\sum_{i=0}^{N-1} x_i)^2}{N(N - 1)}} \\ &= \frac{1}{\sqrt{N(N - 1)}} \sqrt{N \sum_{i=0}^{N-1} x_i^2 - (\sum_{i=0}^{N-1} x_i)^2} \end{aligned} \quad (8)$$

(https://en.wikipedia.org/wiki/Algorithms_for_calculating_variance).

The acquired squared data sq_{bin} returned from the transient recorders is

$$sq_{bin} = \sum_{i=0}^{N-1} x_i^2. \quad (9)$$

With this and the acquired data from eq. 4 the right square root in eq. 8 can be written as

$$sqd_{bin} = \sqrt{N \sum_{i=0}^{N-1} x_i^2 - (\sum_{i=0}^{N-1} x_i)^2} \quad (10)$$

$$= \sqrt{N sq_{bin} - (c_{bin})^2}. \quad (11)$$

This will fit into 4 byte numbers used in the standard data files and will therefore be saved by TCPIP Acquis.

Reading sqd_{bin} back from the data files will enable to restore the *sample standard deviation* s and the *standard error of the mean* σ_{μ} as

$$s = \frac{sqd_{bin}}{\sqrt{N(N-1)}} \text{ and} \tag{12}$$

$$\sigma_{\mu} = \frac{s}{\sqrt{N}}. \tag{13}$$

In TCP/IP Acquis and the Viewer software data according to eq. 13 is displayed.


5.5 Acquisition Low Level Description

Once the Transient Recorder is started it will wait for a trigger pulse on either of its trigger input - A or B. The starting consists of two operations which can also be executed separately: clearing the memory and arming the Transient Recorder. While arming the Transient Recorder is nearly instantaneous, clearing the memory however requires a time similar to that required for an acquisition from memory A and B. For a TR20-160 this would be 6ms (2*3ms).

Once the TR is armed it waits for next trigger. The input where it arrives determines which summation memory will be used. This will continue until either the shotnumber reaches 4094 or a stop command is send. The stop command will not be executed immediately but rather tells the system to return to the idle state and not the armed state after the next acquisition.

Once the system is in the idle state the acquired data can be transferred to the PC. There is no acquisition possible in parallel with data transfer. So this adds to time when the Transient Recorder cannot average. The typical transfer time for DIO-32HS based system is 20ms for a single dataset with 16k. For a Ethernet based system the Ethernet Controller will transmit at 200 Kbytes/second. A single dataset with 16k will then take $32/200 = 160ms$.

Once all the data is transferred the Transient Recorder can be started again.



Acquisition state transition diagram

Shot number considerations The reported shot number is zero after a reset. After the memory has been cleared the TR reports a shot number of two, after the first acquisition a shot number of 3 is reported and so on till 4096 is reported. At this stage the above mentioned 4094 shots are acquired and the TR will stop.

To stop the TR at a predefined shot number below this, one should wait till the desired `shot number + 1` is reached and then issue a stop command. If a trigger is still supplied after this the TR will return from the armed state with the next trigger and the reported shot number will be `shot number + 2`, which consists of the two additional cycles from the memory clear and the desired shot number.

Calculation of lost shots


1. Stratospheric system, lets assume a system with 4 channels TR20-160, both analog and photon counting data needs to be transferred for the whole range(120km). The Laser frequency is 30Hz. The system will need 6ms for initializing all 4 TR, as this can be done in parallel. Then it will acquire 4094 shots and all data is to be transferred to the PC. So one has the LSW and the MSW for the Analog and the PHOTON dataset for the photon counting. That's 60ms per Transient Recorder and 240ms overall for the system. So one would loose 250ms for reading and restarting or 7 to 8 shots per 4094 shots.

For a tropospheric system this would be 96k per Transient Recorder and 400k for the whole system this would require 2sec and one will loose 60 shots.

2. Tropospheric System, assume 1 channel TR20-160, where only the first 1000 bins of analog data are transferred (a 7.5km trace). Then even with the Ethernet Controller the dataset would be transmitted within 10ms to the PC, so no shot would be lost.

Please note that these times assume that the timing constraint comes from the Transient Recorder or the Ethernet Controller. In real systems the PC is also a limiting factor. For instance the start of the data transfer requires to start a DMA process which consumes time down to a couple of ms or the front panel activity like displaying the data may require significant time.

Push Mode For Ethernet based systems sending a lot of small commands can cause delays due to the Nagle Algorithm (see http://en.wikipedia.org/wiki/Nagle's_algorithm for the details). To overcome this the Ethernet Controller is implemented with a "PUSH Mode" mechanism that is activated when for a predefined shotnumber has to be acquired. When this mode is active the Ethernet Controller sends the data (push) without further request to the PC and restarts the Transient Recorder. This will continue till the push mode is revoked. After the push mode is revoked the Controller returns into the normal mode (slave).



Push Mode state transition diagram

5.6 High Performant Read Operations

While the MPUSH mode is the most efficient way to transfer data it creates a burden on the programming as the data is flowing asynchronous to the program flow which creates problems when synchronous action is required like moving a alignment mirror.

The [GetDataSet](#) has the model that it first requests data and then waits till all the data is available. This is again sending small packages and invokes the Nagle algorithm.

By splitting the function in to two functions

- [RequestDataSet](#)
- and using the [Licel_TCPIP_ReadData](#)

The requests can be bunched together which will defeat the Nagle algorithm.

The [ReadSpeedDemo](#) measures the improvement by calling the [GetDataSet](#) function 100 times and then doing 100 requests and 100 retrieve operations.

The speed improvement is heavily dependent from network speed and router transmission and varies from being two times better than the [GetDataSet](#) function to 30% worse.

6 Network security

The Licel Ethernet Controller might be the target of an attack. The best protection against this is to run the Controller with a private IP address beyond a firewall. Firewalls are designed to protect against various types of attacks that can not be covered by the Ethernet Controller. Licel strongly recommends the use of a firewall/router combination to prevent unauthorized use of the hardware.


Starting with firmware versions from 2005-02-22 ([state53](#)) the Licel Ethernet Controller has an additional level of security that can be used additionally.

A secure mode combines white listing of allowed hosts with an encrypted password transmission scheme. In order to activate the secure mode,

- One needs to transmit the white listed hosts to the Controller, and send a connection password when activating the secure mode.
- Once this is done the Controller will check whether the host is authorized to access it.
- It then send a token that the host needs to encrypt with the connection password and send back the encrypted token to the Controller.
- The Controller then decrypts the received encrypted token with the previously received connection password and compares it to what was sent by it for encryption.
- If the response is correct the connection is allowed to proceed otherwise the connection is closed.

The idea behind this is that the connection password is transmitted to the Controller only in a secure environment and later the password is used to encrypt and decrypt a random token.

The algorithm for encryption/decryption is a blowfish algorithm which is a open algorithm without license restrictions. See [Bruce Schneiers Page](#) for the details.



Secure Mode - State transition diagram

The setting of the secure mode will persist during power off and on.

As a default the secure mode is disabled and all hosts can access the Controller. The hardware reset will also reset the secure mode and remove all information about the white listed hosts. If during the secure mode activation something goes wrong, like the controlling PC is not white listed, then the only way to get again access to the Controller is a hardware reset. Due to this it is highly recommended that the secure mode is only enabled when one has physical access to the Controller.

7 C - Example Programs

All example programmes can be compiled under Windows from a developer prompt for Visual Studio. The `makeall.bat` script issues the required `nmake` commands.

Under Linux a `sh buildall.sh` from a terminal prompt will do the same.

This section contains 4 sample program which could be used in a batch file. They demonstrate the basic actions to run the Transient Recorder.

`start.c` - Configures the Transient Recorders for a acquisition and starts them.

`shot.c` - Shows the number of shots already acquired.

`read.out.c` - Transfers the data from the Transient Recorders to the PC and writes them to a data file.

`Show.tmp.c` - Displays the data. The example programs can be compiled by

- `make -f start.mk`
- `make -f shot.mk`
- `make -f read.out.mk`
- `make -f show.tmp.mk`

Under Windows the executable files (.exe) of the examples can be run from the standard Command Line Interface (CLI). For the DIO32HS version adjust the include path and the path to `nidaq32.lib`. Please note that currently there is no routine for `show.tmp` under WinXX.


7.1 The script

```
./start
while true
./shot
./readout
./start
./show_tmp
done;
```


The basic working of the examples provided are as follows :

The [configuration file](#) contains all the Transient Recorder parameter information. This must be configured correctly by the user before executing the programs. The start program reads the configuration file and extracts the necessary parameter information to configure the Transient Recorder(s). The Transient Recorders are then started. Depending on the user configuration of the maximum shots to acquire in the configuration file the 64k shot mode is activated or deactivated. If not activated the Transient Recorder acquires 4094 shots (default value in 4k mode). If activated it acquires 65536 shots (default value in 64k mode). Then the data is transferred and the Transient Recorder(s) is (are) restarted for the next acquisition. While the data is acquiring the previously acquired data is shown to the user. The data is written to a mixed ASCII-binary format data file. The [header](#) of this file contains some key important information about the Transient Recorder and the datasets.


7.2 start



7.3 shot



7.4 readout



7.5 Example Program Configuration File

The examples use a ASCII based configuration file `standard.cfg`. The information is classified into two groups. First the measurement location/situation and second, the configuration info for each dataset.

```
char [8]           Measurement site integer altitude above sea level [m]
double            Longitude
double            Latitude
```

<i>char</i>	<i>Leading letter of filename</i>
<i>char [250]</i>	<i>Output directory for data, must be identical with input directory in mega.cfg</i>
<i>bool</i>	<i>Does the laser shoot while reading, i.e. is there a trigger pulse while reading (no - 0, yes - 1) : range (0-1)</i>
<i>integer</i>	<i>Identification number of the National Instruments board;</i>
<i>integer</i>	<i>Number of data sets</i>
<i>integer</i>	<i>Transient Recorder number (range 0-15)</i>
<i>string</i>	<i>Transient Recorder type in the format TRXX-YY, where XX is the sampling rate. YY can either be the number of ADC Bits of the Transient Recorder ("12bit" or "16bit") or the memory length. The later works for older TR where the default of 12 bit is assumed.</i>
<i>integer</i>	<i>Memory MEMORY_A - 0, MEMORY_B -1 (range 0-1)</i>
<i>integer</i>	<i>Signal type: analog - 0, photon counting - 1 (range 0-1)</i>
<i>integer</i>	<i>Bins number of data bins in the data file, these bins may incorporate more than original Transient Recorder bins, if the data reduction below is larger than 0.</i>
<i>integer</i>	<i>Signal range Analog signal type: 500mV - 0, 100mv - 1, 20mv -2 (range 0-2) Photon counting signal type: discriminator level (range 0-63) 63-1.25V</i>
<i>integer</i>	<i>Show overflow values.</i>
<i>integer</i>	<i>Voltage at photomultiplier [V].</i>
<i>double</i>	<i>Laser frequency [Hz].</i>
<i>double</i>	<i>Number of bins.</i>
<i>integer</i>	<i>Data reduction factor (2^n), this for data reduction of 2 for instance $2^2 = 4$ Transient Recorder bins will be combined into a single data bin.</i>
<i>integer</i>	<i>Polarization none -0, parallel -1, perpendicular -2 double wavelength [nm]</i>
<i>integer</i>	<i>Laser source identifier</i>

7.6 File format

The example program uses the same file form at as the LabVIEW `TCPIP-Acquis.llb`. By this method the files are inter operable between the different platforms. The file format is a mixed ASCII-binary format where the first lines describe the measurement location/situation, below follow the dataset description and then finally the raw data as 32-bit integer values itself.

Sample file header

```
a9981017.204567
Berlin 10/08/1999 17:20:36 10/08/1999 17:20:41 0015 0015.0 0053.0 00
0000000 0010 0002000 0005 02
1 0 2 08000 1 1600 07.5 286.0 0 0 00 000 12 002000 0.100 BT1
1 1 2 08000 1 1600 07.5 286.0 0 0 00 000 00 002000 0.793 BC1
```

Line 1

Filename *string a.*
Format: ?yyMddhh.mmsmsms

? - The first letter can be choose freely.
yy - Two numbers showing the years in the century.
M - One number containing the month as a hexadecimal number
(December ≡ C).
dd - Two numbers containing the day of month.
hh - Two numbers containing the hours since midnight.
mm - Two numbers containing the minutes.
s - Two number containing the seconds.
ms - Two number containing the milliseconds divided by ten.

Line 2

Location *String with 8 Letters.*
Start Time *dd/mm/yyyy hh:mm:ss..*
Stop Time *dd/mm/yyyy hh:mm:ss.*
Hight asl. *Four digits (meter).*
Longitude *Four digits (including - sign). one digit for decimal grades.*
Latitude *Four digits (including - sign). one digit for decimal grades.*
zenith angle *Two digits in degrees.*

Line 3

Laser 1 Number of shots *Integer 7 digits*
Pulse repetition frequency for Laser 1 *Integer 5 digits*
Laser 2 Number of shots *Integer 7 digits*
Pulse repetition frequency for Laser 2 *Integer 5 digits*
number of datasets in the file *Integer 2 digits*

Dataset Description

Active *1 if dataset is present, 0 otherwise*
Analog/Photoncounting *Analog ≡ 0, Photoncounting ≡ 1*
Laser source *One digit Laser 1 ≡ 1, Laser 2 ≡ 2.*
Number of bins *5 digits*
1
PMT highvoltage *Four digits in Volt*
binwidth *In meter two digits before . and 2 digits after the dot*
Laser wavelength *In nm, three digits dot*
Polarisation *One letter, o ≡ no polarisation, s ≡ perpendikular, l ≡ parallel*

<i>0 0 00 000</i>	<i>Backward compatibility</i>
<i>number of ADC bits</i>	<i>In case of an analog dataset, otherwise 0</i>
<i>number of shots</i>	<i>6 digits</i>
<i>analog input range/discriminator level</i>	<i>Analog input range in Volt in case of analog dataset , discriminator level in case of photon counting, one digit dot 3 digits.</i>
<i>Dataset descriptor</i>	<i>BT ≡ analog dataset, BC ≡ photon counting, the number is the transient recorder number as a hexadecimal.</i>

The data set description is followed by an extra CRLF. The datasets are 32bit integer values and are separated by CRLF. The last dataset is followed by a CRLF. These CRLF are used as markers and can be used as check points for file integrity.

7.7 SampleAcquis for Ethernet Applications


This is a skeleton for a simple acquisition with a Ethernet system.



7.8 MPUSH Mode Sample Acquisition


This is a demo how to use the MPUSH Mode to store data. By giving the Ethernet controller the control of the acquisition cycle the data throughput is optimized. The MPUSH mode gives the full potential if data sets with a

small number of shots are frequently transferred. All the functions from start, shot, and readout are combined here.




7.9 APD Control example

In the `SetAPDHV_src` is the code `setAPDHV.cpp` which reads the command line and sets the TEC controller state and waits till the TEC is stabilized and sets then the HV.



7.10 Wind Sample Acquisition

This is a demo how to use the Wind Acquisition to store data. All the functions from `start`, `runs`, `shot`, and `readout` are combined here.



7.11 Network management utilities

7.11.1 Getting Started

This module shows the basic process of connecting to a Licel Ethernet Controller. The Controller will return the ID and reveal the capabilities that can be used with this Controller.

7.11.2 Set Fixed IP Address

This module sets the Controller IP address. The new address will be activated after the Controller is turned off and on again

7.11.3 Activate DHCP Mode

This module activates the DHCP mode. The Controller will enter DHCP mode once it is turned off and on again

7.11.4 SecureModeEnable

This module enables the secure mode for accessing the Ethernet Controller. See the [Network Security](#) section for more details

7.11.5 SecureModeDisable

This module disables the secure mode for accessing the Ethernet Controller. The Controller will be fully accessible for all hosts. See the [Network Security](#) section for more details

8 Appendix VB6/VB.net Programming

The Visual Basic driver is not part of our standard distribution and can be ordered from Licel separately. Licel recommends the use of the LabVIEW Modules, which allows a much faster programming and deployment of your applications. The VB.net driver should be used when existing applications are upgraded to communicate with the Licel Ethernet Controller. The VB6 programming is provided for legacy programming when the upgrade to VB.Net is not possible. The programs require the Installation of SocketWrench. The installer is provided as a zip-file in the VB6 directory. The VB6 and the VB.net use a control that ships with LabVIEW to display the data. If you do not LabVIEW you will need to replace the control first before running the applications.

8.1 Sample applications

For demonstration purposes some demo modules are available. The main purpose is to show the use of the functions provided in the driver, which are described after the modules.

8.1.1 Control Overview

This module demonstrates the use of the different commands for controlling a single Transient Recorder, PMT's, APD's and the trigger module. When the module is started and the connection is established the capabilities of the Controller are queried. For each supported capability a tab will become available

8.1.2 MultipleChannel

This module demonstrates the use of commands to control multiple transient recorder with single commands. This might be necessary for performance reasons as sending a lot of small commands and waiting every time for the response can yield large delays as the OS tries to delay the sending of TCPIP packages that are almost empty. The starting point is the `Start_click` function. It opens a connection to the Controller and starts the selected TR's then it launches a timer which will every second readout the TR's. If the selected TR's have reached 4094 shots the memory will be cleared and the acquisition will be restarted.

8.1.3 PushModeDemo

This module demonstrates the use of the push mode for fast data transfer. The push mode can be used when repetitively data needs to be transferred between the TR and the PC. If the shot number is low this would require a lot of calls to StartAcquisition, GetStatus and GetDataSet. In the push mode one sends only one command at the start (see the Start_Click routine). Then a timer is started and the data is read as it arrives. Please note that the data is transferred over a second socket connection.

9 Appendix - Obsolete DIO32HS references

This is held only for reference if you encounter an old DIO based system. The way to handle this is to upgrade the system to an Ethernet based system, which can be done also for very old systems.

9.1 Operation Principles

The communication with the Transient Recorder is performed via a parallel bus using a hardware handshake. The parallel bus is based on the digital I/O card family DIO-32HS supplied by National Instruments. For the description of the bus timing refer to the [user manual](#) from National Instruments. The used protocol is Level-Acq. This interface is also implemented internally between the Ethernet Controller and the Transient Recorders.

9.2 Hardware Requirements

Operating environment	Card type	Required Slot
PC	Ethernet deprecated	Ethernet connector
Desktop-PC Win9xx	AT-DIO-32HS	ISA Slot short
Win-NT	PCI-DIO-32HS	PCI slot
Notebook	DAQCard 6533	PCMCIA slot
PXI	PXI 6533	PXI slot

9.2.1 Further References

1. [DIO-653x User Manual](#)
2. [NI-DAQ User Manual for PC Compatibles](#)
3. [NI-DAQ Function Reference Manual for PC Compatibles](#)
4. [LabVIEW Measurements Manual](#)

9.3 Windows - DIO-32HS

9.3.1 Required software

New systems should not be installed with this option. This will go into the unsupported state as Windows XP market share becomes negligible.

In order to use the C, Basic or LabVIEW routines under windows you will need a working installation of

- Visual C++ 6.0 or
- Visual Basic 6.0 .
- LabVIEW 7.x or higher

9.3.2 NI-DAQ-Setup

- Check first whether you find `Measurement and Automation Explorer` in your National Instruments program group or a `NIDAQ configuration Utility` in your LabVIEW program group.
- If present start them and look for the NIDAQ-Version. We found that the most reliable version are 5.1.1, 6.5.1 and 6.9.3
- The NIDAQ-Versions 6.0 - 6.1 do not work properly with the interface card, if you have these versions already on your PC make sure to change to 5.1.1 or a higher version
- There have been reports about crashes with NI-DAQ 6.8.x under Win98.
- If not present install the NIDAQ-Software, preferable 6.9.3, as indicated in the NIDAQ Manual.
- In order to compile the corresponding projects the language support for Visual C or Visual Basic should be installed.
- Once the language support is installed there should be a `nidaq.h` in `Program\National Instruments\Ni-daq\Include\` and a `nidaq32.lib` in `Program\National Instruments\Ni-daq\Include\` for a working C installation.
- For Visual Basic you should be able to locate `nidaq32.bas` in the `Program\National Instruments\Ni-daq\Include\` directory.

9.3.3 Interface card installation

The DIO-32HS cards are Plug and Play enabled cards. Windows9x detects the plug and play cards during the boot process. There are special considerations with ISA PNP boards under Windows NT.

Configuring ISA Plug and Play Devices for Windows NT 4.0 If you plan to use ISA Plug and Play DAQ devices on Windows NT 4.0, you must first install the Windows NT 4.0 ISA Plug and Play driver before configuring your device with the NI-DAQ Configuration Utility. This driver is not installed by default. Follow these steps to install the driver:

1. Insert your Windows NT 4.0 CD.
2. Go the `\Drvlib\Pnpisa\X86` directory.
3. Right-click once on the `Pnpisa.inf` file, select the Install option, and follow the instructions.
4. After you have installed the `Pnpisa.inf` file, shut down your computer.
5. Install your ISA Plug and Play DAQ device.
6. Turn on your computer. When Windows NT 4.0 detects your ISA Plug and Play DAQ device, it will specify the necessary driver files. Because this will result in a configuration change, restart your computer.
7. After you have restarted your computer, run the NI-DAQ Configuration Utility to configure your device.

Configuring PCI Boards

1. PCI boards will be detected during the boot process.
2. Turn the computer off
3. Insert the card and boot the machine
4. The card will be detected during the boot process.
5. The following resources are necessary
 - I/O Ports
 - one Interrupt
 - one DMA Channel
6. The configuration can be changed under the Windows9x in the system manager and under Windows NT with the NIDAQ Configuration utility.

9.3.4 Verification

Assuming that you have managed to bring up the machine again without any dirty messages about resource conflicts, you have to test whether the NIDAQ software recognizes the DIO-32HS.

Remove all cables from the DIO-32HS.

Open the Measurement and Automation Explorer. Press the Testpanel button. If you have resource conflicts you should go back to the hardware manager and change the configuration there, reboot and repeat the verification. If your configuration is valid, select the same port for input and output. After setting some lines in the output to TRUE you should see the same lines true in the input port, because the card reads back the output lines. If it does not work you have to change the configuration again and repeat the procedure above.

9.4 Linux - PCI-DIO-32HS

*New systems should **NOT** be installed with this option. This is now in the unsupported state as the hassle of kernel compilation does not give any advantage over the Ethernet option. Older systems can be upgraded to the Ethernet option. The following is documented for archive references purpose only.*

9.4.1 Kernel preparation

The COMEDI driver is compiled versus the kernel sources. In order to load the modules correctly, the kernel sources should correspond to the booted kernel. The best way to ensure this is to compile a new kernel and make it bootable. Instruction on doing this are beyond the scope of this manual. For details please refer either the [Kernel HOWTO](#) or the vendor documentation of your distributor. Please note that every time you update the kernel you will need to recompile comedi, comedilib and all other applications.

9.4.2 Necessary Files

The modified COMEDI versions are distributed as two archives comedi.tgz and comedilib.tgz. Please create first a directory, for example `com` and copy both the archives there. Unpack them by

- `tar -xzvf comedi.tgz`
- `tar -xzvf comedilib.tgz`

The following files have been modified with respect to the original version.

- `comedi/comedi/comedi_fops.c`
- `comedi/comedi/kvmmem.h`
- `comedi/comedi/drivers/ni_pcidio.c`
- `comedi/comedi/include/linux/comedi.h`
- `comedilib/include/comedi.h`
- `comedilib/include/comedilib.h`
- `comedilib/lib/Makefile`
- `comedilib/lib/dio.c`
- `comedilib/lib/dio_licel.c`

9.4.3 Driver tests and card installation

Please follow the instructions outlined at `comedi/INSTALL` and `comedilib/INSTALL` to install both packages.

Once both packages are installed and `/etc/modules.conf` modified, one should be able to run `/sbin/modprobe ni_pcidio` without any errors.

Shutdown the computer, plugin the PCI-DIO-32HS and reboot the machine.

Run again `/sbin/modprobe ni_pcidio` and `/usr/sbin/comedi_config /dev/comedi0 ni_pcidio`. It should not generate an error. Make sure to check `/var/log/messages`

9.4.4 Changes to /etc/modules.conf

The following to lines should be added to /etc/modules.conf:

```
alias char-major-98 comedi
alias char-major-98-0 ni_pcidi0
```

9.4.5 Changes to /etc/rc.d/rc.local

The following to lines start the module at every bootup. Otherwise these commands should be issued by the superuser:

```
/sbin/modprobe ni_pcidi0
/usr/sbin/comedi_config /dev/comedi0 ni_pcidi0
```

9.4.6 Directory structure

The comedi driver is placed below the com2 directory. The C-Sources for the examples are below licellinux. The LabView Libraries are in the licellinux/labview/. The code for the interface between LabView and comedi is in the licellinux/labview/lib folder.

9.4.7 LabView + comedi

Recompiling comedilib may require a recompilation of the glue code. This can be done by

- make clean
- make

In the licellinux/labview/lib/config, ..read and ..write directories. Once the CIN code resources are rebuild they should be reloaded into the corresponding CIN's in licellinux/labview/LV.COMEDI.Interface.llb. In the following VI's a code resource reload would be necessary:

- lv_config.vi
- lv_read.vi
- lv_write.vi

9.4.8 licel_nidaq/licel_re

The licel_nidaq.(c|bas)/licel_re.c are used for the National Instruments interface cards and have four routines inside:

```
/* Setup of the DIO card*/
int InitializeBoard(short int iBoard)
Public Function InitializeBoard(iBoard As Integer) As Integer

/* Writing the commands to the TR*/
int WriteCommand(short int iBoard, short iCommand)
Public Function WriteCommand(iBoard As Integer, iCommand As Integer)
As Integer

/* Reading back information from the TR*/
int ReadArray( short int iBoard, unsigned short *piBuffer,
unsigned long ulCount)
Public Function ReadArray(iBoard As Integer, piBuffer() As Integer,
ulCount As Long) As Integer

/* Selection of TR group*/
int Select(short int iBoard, int hilow)
Public Function Select.TR(iBoard As Integer, iDevice As Integer) As Integer
```

In general these routines do not need a customer tweaking, the only place where customization may be useful is in InitializeBoard the `AcqDelay` can be changed between 0 (0ns) and 7 (700ns), which might help if data transfer errors do occur (readout errors usually give variable readout while reading repetitive the same dataset, if this happens it is the right time to contact Licel for help)